

NeuroPapain™ Enzyme

Neuronal Tissue Dissociation Enzyme



A division of Gene Therapy Systems, Inc.

Contents		
Cat. No.	Description	Quantity
NM100200	NeuroPapain Enzyme	100 mg

Related Products	Cat. No.
NeuroPrep™ Medium, 100 ml	NM100100
NeuroPure™ E18 Primary Rat Hippocampal Cells	N100200
NeuroPure™ E18 Primary Rat Cortical Cells	N200200
NeuroPure™ P8 Primary Rat Cerebellar Cells	N300200
NeuroPure™ Primary E18 Hypothalamus Cells	1 pair
NeuroPure™ Primary E18 Striatum Cells	~1 x 10 ⁶ cells*
NeuroPure™ Primary E18 Spinal Cord Cells	~1 x 10 ⁶ cells*
NeuroPure™ Primary E18 Midbrain Cells	1 pair

Shipping & Storage	NeuroPapain Enzyme is shipped at room temperature and should be stored at 4° C upon receipt. It is stable for 6 months when stored properly.
-------------------------------	--

INTRODUCTION

The NeuroPapain™ Enzyme is an optimized formulation of the enzyme papain that is ideal for enzymatic dissociation of neuronal tissues such as the NeuroPure and NeuroStem™ Primary Rat Neurons. NeuroPapain Enzyme provides rapid and safe enzymatic digestion of primary rodent neuronal tissues when used in conjunction with NeuroPrep Medium (Cat#: NM100100). Following treatment with NeuroPapain in NeuroPrep Medium, neuronal tissues are more easily dissociated via mechanical trituration, and the number of viable dissociated cells can be increased by up to 100%. However, please note that for assays performed within 4 days of plating, some digestion of surface proteins is inevitable.

MATERIALS AND METHODS

Note: The following protocol has been validated for use with the NeuroPure Primary Rat Neuronal Cells. Please refer to the NeuroPure protocols for complete details.

1. Add 5 mg of NeuroPapain Enzyme into 2.5 ml of NeuroPrep Medium. Mix at 37°C for 15 minutes to completely dissolve the NeuroPapain. Sterilize this solution with a 0.2 µm filter prior to utilizing for tissue digestion. Use within 3 hours for best results.
2. Prior to enzymatic treatment, allow the neuronal tissue to settle for 15 - 30 minutes at 4°C. Alternatively, place the tube containing the tissue in a 50 ml tube and spin down the cells at 1,100 rpm (200xg) for 1 minute. Transfer the medium from the cell vial to a separate sterile tube while being careful not to remove any loose tissue pieces. Save the medium for trituration following NeuroPapain treatment.
3. Immediately add 2 ml of sterile NeuroPapain solution to the tissue-containing tube, and allow the neuronal tissue to incubate for 30 minutes at 30 °C. Swirl every two minutes by hand.
4. Following incubation, spin down the cells at 1,100 rpm (200xg) for 1 minute. Remove the NeuroPapain solution, again being careful not to disturb or remove the tissue.
5. Add 1 ml of shipping medium back to the neuronal cells. Save the other 1 ml of shipping medium for Step 4 in the NeuroPure protocol.
6. Proceed to Step 3 in the NeuroPure protocol.

Mike Vengrow 12/5/06 10:24 AM

Deleted: [NOTE: you don't need to use the "TM" mark each time you mention the product name; the convention is to use the TM indication only the first time you mention the product on each page.]