

Transformation in 3 Minutes: TurboCells™ Chemically Competent *E. coli*

TurboCell™ Chemically Competent *E. coli*

Questions & Answers

1. What is the primary difference between the TurboCells™ Chemically Competent *E. coli* and standard competent *E. coli*?

The TurboCells Chemically Competent *E. coli* are specially optimized to achieve high transformation efficiency with a novel 3-minute transformation protocol. In contrast, transformation using standard competent cells typically requires a 1.5 – 2 hour protocol. Figure 1 demonstrates the efficiency of the TurboCells protocol compared to a standard transformation protocol.

2. What transformation efficiencies can be achieved using TurboCells Chemically Competent *E. coli*?

When using the 3-minute TurboCells transformation protocol, you can achieve efficiencies of 5×10^7 - 1×10^8 cfu/ μ g of supercoiled DNA. Also, when using a traditional 1.5 - 2 hour protocol with the TurboCells, greater than 1×10^9 cfu/ μ g of supercoiled DNA can be achieved.

3. What is the difference between TurboCells™ and TurboCells F'?

TurboCells F' contain a self-transmissible, low-copy plasmid, called F', that can be used for the generation of single strand DNA (e.g., using M13 bacteriophage techniques).

4. In addition to time-savings, what other advantages do TurboCells Chemically Competent *E. coli* offer?

TurboCells Chemically Competent *E. coli* are prepared according to a unique procedure that allows robust performance under normally deleterious transformation conditions. For example, it is not required to dilute or purify your ligation mix before transformation. If needed, over 10 μ l of full strength ligation mix can be added to 50 μ l of competent cells without significantly compromising transformation results. The genotype of TurboCells is suitable for most cloning needs, such as blue/white screening, generation of plasmid based libraries, and the ability to efficiently transform large plasmids.

5. What are the genotypes of the TurboCells Chemically Competent *E. coli*?

The genotypes of the TurboCells are summarized in the table 1.

Figure 1: Comparison of TurboCells™ Chemically Competent *E. coli* Protocol and Standard Chemically Competent *E. coli* Protocol




Figure 2. Transformation Results with TurboCells



pUC19 + TurboCells
on IPTG/Amp plates

pUC19 self ligation+
TurboCells on
IPTG/Amp plates

6. Why are TurboCells Competent *E. coli* packaged in single-use aliquots?

TurboCells are packaged in single-use aliquots to avoid efficiency-robbing freeze/thaws and wasted cells.

7. Are TurboCells available for protein expression?

Yes. TurboCells BL21(DE3) and BL21(DE3) pLysS offers 3-minute transformation of T7-promoter based expression vectors for protein expression.

8. What concentration of kanamycin and ampicillin should I use with TurboCells?

We use 50 µg/ml and 100 µg/ml, respectively.

Table 1. Genotype of TurboCells™ Chemically Competent *E. coli*

TurboCells™: *recA1 endA1 hsdR17 supE44 thi-1 gyrA96 relA1 φ80lacZΔM15 Δ(lacZYA-argF)U169*

TurboCells™ F': *F recA1 endA1 hsdR17 supE44 thi-1 gyrA96 relA1 φ80lacZΔM15 Δ(lacZYA-argF)U169*

TurboCells™ BL21(DE3): *recA1 endA1 hsdR17 supE44 thi-1 gyrA96 relA1 φ80lacZΔM15 Δ(lacZYA-argF)U169 (DE3)*

TurboCells™ BL21(DE3)pLysS: *recA1 endA1 hsdR17 supE44 thi-1 gyrA96 relA1 φ80lacZΔM15 Δ(lacZYA-argF)U169 (DE3)pLysS(Cam^R)*

| Genotype | Advantage |
|-----------------|--|
| <i>recA1</i> | Mutation in gene(s) responsible for recombination of DNA. This genotype is particularly desirable when cloning genes with direct repeats. |
| <i>endA</i> | Mutation in the nonspecific endonuclease Endonuclease I. Eliminates non-specific endonuclease activity resulting in improved plasmid preps. |
| <i>hsd</i> | Mutations in the system of methylation and restriction which allows <i>E. coli</i> to recognize DNA as foreign. The <i>hsd</i> genotype allows efficient transformation of DNA generated from PCR reactions. |
| <i>lacZΔM15</i> | Element required for β-galactosidase complementation when plated on X-Gal. Used in blue/white screening of recombinants. Usually carried on the lambdoid prophage 80 or F'. |
| <i>DE3</i> | Lysogen that encodes T7 RNA polymerase. Used to induce expression in T7-driven expression systems. |
| <i>pLysS</i> | Plasmid that encodes T7 lysozyme. Used to reduce basal expression in T7-driven expression systems by inhibiting basal levels of T7 RNA polymerase. |

| Product | Quantity | Catalog no. |
|---|------------|-------------|
| TurboCells™ Chemically Competent <i>E. coli</i> | 20 x 50 µl | C300020 |
| TurboCells™ F' Chemically Competent <i>E. coli</i> | 20 x 50 µl | C301020 |
| TurboCells™ BL21 (D3E) Chemically Competent <i>E. coli</i> | 20 x 50 µl | C302020 |
| TurboCells™ BL21 (D3E) pLysS Chemically Competent <i>E. coli</i> | 20 x 50 µl | C303020 |