STCON3 Conductivity Measuring Cell
Celda de medición de conductividad STCON3
Cellule de mesure de la conductivité STCON3

Ohaus Corporation
7 Campus Drive
Suite 310
Parsippany, NJ 07054 USA
Tel: (973) 377-9000
Fax: (973) 944-7177

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This Conductivity Measuring Cell utilizes the 4-Ring potentiometric method for measuring conductivity, which incorporates a series of four stainless steel rings formed into the probe shaft.

This design completely eliminates the polarization, which occurs with the 2-Plates amperometric method.

Furthermore, without polarization the probe can measure a wider range of conductivity values because it does not suffer from electrolysis.

Order Numbers:
83033972

Operating Instructions

1. Make sure the plastic shield is in place when measuring.
2. When measuring make the solution reach the line on the plastic shield and below the vent hole.
3. To prevent carry over from high to low conductivity solutions, rinse with distilled water between measurements.
4. Make sure the cell chamber is bubble-free when measuring. To reduce air bubbles, immerse probe in the solution at an angle and then rise to a vertical position.
5. Allow sufficient time for the sensor to stabilize when measuring samples at different temperatures. Manual end-pointing is advised.

6. The sensor is not recommended for low ionic strength solutions. (<70 μ S/cm)
7. Clean the probe and shield with distilled water after use.
8. For Calibration and measurement procedures, refer to Conductivity Meter Instruction Manual.

Reagents and Supplies :

1413 μ S/cm Conductivity Standard
12.88mS/cm Conductivity Standard

Precautions and Limitations :

1. Always have the plastic shield in place when measuring.
2. Don't expose the probe to organic solvents when either cleaning or taking measurements.
3. Don't use the probe outside the recommended temperature range of 0 to 50°C.
4. Specification and performance of this probe is guaranteed only when used with the appropriate model of Ohaus Conductivity Meter.

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Ohaus Corporation, all rights reserved/ todos los derechos reservados/ tous droits réservés