General instructions for STREF electrodes

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0.198 V (STREF1)

0.241 V (STREF2)

A reference electrode is an electrode which has a stable and well-known electrode potential.

The electrode is delivered ready for use. Please check visually for broken parts.

Connect the reference electrode to the reference input(2mm banana) socket of your pH meter(ST2100 or ST3100). Connect the measuring electrode to the appropriate input of your meter.

Clutch the reference electrode and shake up and down in order reduce bubbles. The electrode is ready for use in conjunction with a working electrode. Please follow the detailed calibration procedure described in the manual of your pH meter.

Storage

Rinse or clean electrode after use and replace the wetting cap, partly filled (1/3) with reference electrolyte.

Refilling:

Refill the electrode with recommended reference electrolyte, if necessary. The electrode is approximately 10mm below the filling hole.

Trouble-shooting

Slow response drift: Check connections; clean diaphragm; after dry storage soak for 24 hours in reference electrolyte.

Slope not adjustable: Clean diaphragm; clean connector

Zero not adjustable: Clean diaphragm; change reference electrolyte.

NOTE:

1. Make sure that the reference electrode compartments are filled with electrolytes solution. (3M KCl)
2. Make sure the junction is not blocked.

The useful life of the reference electrode can be maximized if used and maintenance properly.

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Ohaus Corporation, all rights reserved/ todos