

e-pill® HIProtector® Voice Medication Reminder and Alert Wellness Wizard (960712)

Overview

Your new e-pill® HIProtector® Voice Medication Reminder and Alert Wellness Wizard (www.epill.com/voicemail.html) is a simple to use but innovative Electronic Caregiver's Assistant. This medication reminder device was developed to enhance the health, safety and wellbeing of people who spend time alone at home (Age-In-Place).



Benefits to the patient and their families / caregivers:

- **Medication Reminder:** Never forget your medications again. Plays customized reminder messages in the voice of a loved one to remind users to take their medications and/or perform other important tasks.

The Wellness Wizard can be programmed to play reminder messages at specified times and dates (daily, every-other-day, weekly, monthly, or any date in the future). In addition to medication reminder messages, other important messages might be: reminders to the user to lock their doors before bedtime, upcoming doctor's appointments, bills to be paid, relatives to call, upcoming birthdays.

- **Activity Monitor:** Indirectly checks on the user's well being.



If patient acknowledges their reminder message, this is an indication that no assistance is needed.

- **Emergency Dialer:** If the user does not acknowledge the reminder message, the Wellness Wizard calls up to 6 designated telephone numbers with an automated alert to inform caregivers that the user may need assistance.

Thus, if the caregivers do not receive an alert call, it is likely that their loved-one is safe and sound.

Quick Set-Up Guide

Installation is as easy as connecting an answering machine. Simply connect the device to a telephone wall jack and to the user's telephone. Then connect the power pack to an electrical outlet. That is it!



The caregiver can easily record / change customized reminder messages in their own voice onto the Wellness Wizard from the telephone connected to the Wellness Wizard *or from any remote telephone* (caregiver's regular home phone or cell phone). Simply listen to and respond to voice prompts.

How It Works

The Wellness Wizard can play up to 20 messages at specified times and dates. When the message plays, the big yellow button blinks. The patient must acknowledge that they have heard the reminder message by pressing this button. If the button is not pressed, the Wellness Wizard repeats the message for 20 minutes and the button continues to blink. If the message is still not acknowledged, the Wellness Wizard makes an automated alert call to the designated caregivers phone or cell phone.



Additional Features

The Wellness Wizard also performs other safety and wellness-related tasks:

- Can be programmed to call designated telephone numbers with reminder messages (example “call CVS pharmacy to reorder Rx prescriptions”).
- Speed dial (4 speed dial buttons) allowing the user to easily call: Caregiver, Doctor, Pharmacy and 911.
- If the user is without power or has used the Wellness Wizard to call 911, the Wellness Wizard will alert the caregivers.
- Voice Clock (by pressing the orange bar, the Wellness Wizard actually tells the user the time, day and date).

**e-pill® HIProtector®
Voice Medication Reminder and Alert Wellness Wizard (Item# 960712)**

Manufactured by:
Safety Response Systems
12 Lake Drive West, Wayne, NJ 07474, USA

Distributed by:
e-pill, LLC d.b.a. e-pill Medication Reminders
70 Walnut Street, Wellesley, MA 02481, USA
www.epill.com or 1-800-549-0095