



BMT-20-S-M-gamma Tube for grinding with stainless steel balls (20x)

Components: Tube, stainless steel balls, cover

Volume: 2 - 15 ml

Application: Dry milling of dry and brittle samples (e.g. kaolin, gypsum, coloured pigments, tablets), cell maceration, processing of materials mixed with fluids

Further suitable for:

- additional media quantities may be added during operation
- sampling during operation
- sample preparation under sterile conditions
- germ free conservation of samples (tissue, blood, etc.)

Packing unit: 20 pieces