

Cryosurgery

Nitrospray Plus

Trigger Operated, Closed System,
Liquid Nitrogen Cryosurgical Unit

- Improved Seal⁶
Eliminates o-ring positioning.
- Grooved Head
Easy to use, comfortable to hold.
- Curved Trigger
Unobstructed view of treatment site.
- Notched Head Grip
Easy to remove/install cap.
- Three (3) year warranty

*30-Day
Unconditional
Warranty*

Static holding time:
8 - 12 hours


Intermittent use:
3 - 6 hours of patient treatment time



6060 Nitrospray Plus
(16oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



6065 Nitrospray Plus Lite
(10oz fill capacity)
Includes 5 spray tips,
tip protector and
base holder



6. Pressure seal ring retrofit kit (1006066) to eliminate o-ring on units prior to RP-11061 and LP-10952 is now available. Not to be used with RP-10454 to RP-10561. Contact Customer Service for assistance.